

## How to Update Registration when you are using Manual Methods

When testing samples using Kits or where the test result doesn't need an instrument to interpret results, you have to register the Manual Determination as your instrument.

### 3-step Process

1. Select the analytes for testing by clicking the box next to it and then click on "Choose Instrument for Selected Analytes".

The screenshot shows a software interface with a top navigation bar containing 'Exit' and 'Sample Conditions' buttons. Below this is a progress indicator with three steps: '1 Registration', '2 Sample Conditions', and '3 Results'. The main content area is titled 'Analytes Requiring Action [1 Analyte]' and contains two buttons: 'Choose Instrument for Selected Analytes' and 'Do Not Test Selected Analytes'. A table below lists analytes with checkboxes:

#	<input checked="" type="checkbox"/>	Analyte
1	<input checked="" type="checkbox"/>	hCG

At the bottom of the interface, there are two buttons: 'Choose Instrument for Selected Analytes' (circled in red) and 'Do Not Test Selected Analytes', and an 'Open Next Instrument' button.

2. Select the Not Applicable as the Manufacturer and Manual Determination as the Instrument Model\* from the dropdown list. Then click on the Add & Select button. This action will take you back to the Registration page.

Choose Instrument for Selected Analytes ✕

Register New Instrument

Manufacturer*	Instrument Model*	Nickname	Serial Number	Action
Not Applicable	Manual determination			Add & Select

Instruments Currently Registered for This Subscription

Instrument Model*	Nickname	Serial Number	Action

Other Instruments

Instrument Model*	Nickname	Serial Number	Action

Close

- 1 Registration
- 2 Sample Conditions
- 3 Results

▼ No Instrument [1 Analyte] ?

▼ Analytes Not Being Tested [13 Analytes] ?

▲ Manual determination [1 Analyte] ?

Change Instrument for Selected Analytes    Do Not Test Selected Analytes

#	<input checked="" type="checkbox"/>	Analyte	Reagent Manufacturer	Reagent Catalog Number	Reagent Submethod
1	<input checked="" type="checkbox"/>	hCG	--Select Manufacturer--	--Select Catalog Number--	--Select Submethod--

Change Instrument for Selected Analytes    Do Not Test Selected Analytes

- Choose from the dropdown list and select. Both Catalog and Submethod information can be found on the Package Insert. If not found on the list, please select the Not listed option to send a request to add it.