How to Change Reporting Units

- a. Go to www.oneworldaccuracy.com and login to your OASYS account.
- b. On the home screen click Test Event Dashboard or Results.

PROFILE	RESULTS	LAB REPORTS	ADMIN REPORTS	NETWORK REPORTS	HELP
What would you lil	ke to do?				
(Test Event Dashboard		Performance Reports	Support Center	
Vhat's new?					
OASYS v15.1 ha	is been released! Click he	re for the Release Notes ar	nd additional information.		

c. Under the current test event click under results next to the program you are changing the units for.

	eworld curacy" Canada					- Home -	Log Out
Test Event Dash	board						
OASYS ID							
				Report D	Delivery Preferences	Test Event Calenda	w - 2015
Test Event Op	en: Jun 17 - Results Deadline: Jun 24 Status: Open						
Subscription	Program	Instructions	Worksheets	Registration	Results	Con	firmation
52545	FOBT4123 Fecal Occult Blood	8	=				8
Test Event Op	en: Sep 30 - Results Deadline: Oct 7 Status: Scheduled						
Subscription	Program	Instructions	Worksheets	Registration	Results	Con	firmation
38707	GFRM733 eGFR Monitoring	-	-	>			

- d. Select the instrument panel that contains the analyte that you would like to change the reporting units for.
- e. Next to the analyte is the column that contains the reporting unit. Use the drop down box to select your reporting unit.

Т	Regi	etration 2 Bample Conc	Beis 🔂 Pa	Agroinstath 🚺 Nooto					U M
A	bott	Cell Dyn 3700 [1 Analyte]	9						
G	alte	LHS00 (4 Analytes) 🛛 🛛							
a	ange	instrument für Selected Analytes	Do Nut Test Ber	eded Analytes					
	11	Analyle	Unit	leagent Blanufacturer		Reagent Catalog Number		Respect Submethod	
	-	CD19 (B Cells)	% positive c +	Beckman Coulter (Coulter Instructs)	•	Coulter EPICS CAT#INIt applicable	•	Opticitation	
	11	CD3 (T Cells)	% positive c +	Beckman Coulter (Coulter Products)	٠	Coulter EPICS CAT# Not applicable		Opticitation	•
	8	CD4 (T Helper)	% positive c +	Beckman Coulter (Coulter Products)		Coulter EPICS CAT# Not applicable		Opticitatier	
	11	CD45+ (Leukocytes) - Counts	6L +	Beckman Coulter (Coulter Products)		Coulter EPICS CAT# Not applicable		OpticLaser	
			4			H.			,
	únge	instrument for Delected Analytes	Do Not Test Set	ected Analytes					
0	en hi	of instrument							

Once the change has been made it will be automatically updated in the system.